

MATERIAL SAFETY DATA SHEET

PROFESSIONAL DISPOSABLES INC.
THE HEALTHCARE DIVISION OF
NICE-PAK PRODUCTS, INC.
TWO NICE PAK PARK
ORANGEBURG, NY 10962-1376
845-365-1700



REVISION DATE: 6/04

I. PRODUCT IDENTIFICATION PRODUCT/TRADE NAME

SUPER SANI CLOTH GERMICIDAL DISPOSABLE WIPE - Q55172, Q43472, H04082, U87295, Q86984

HAZARD RATING(NFPA)

HEALTH: 1
FLAMMABILITY: 3
REACTIVITY: 0

SPECIFIC: NONE

EMERGENCY OR INFORMATION TELEPHONE NO: 845-365-1700 (M-F DAYTIME) AT OTHER TIMES, CONTACT THE LOCAL POISON CONTROL CENTER

CHEMICAL NAME: PROPRIETARY MIXTURE

II. HAZARDOUS INGREDIENTS

PER 29 CFR 1910.1200

Hazardous Ingredients	%	ACGIH TLV	CAS NUMBER
ISOPROPYL ALCOHOL	55	400 ppm	67-63-0
QUATERNARY AMMONIUM CMPD	<1	N/A	68391-01-5/68956-79-6

III. PHYSICAL/CHEMICAL CHARACTERISTICS

COLOR/ODOR/APPEARANCE: APPLICATOR SATURATED WITH CLEAR LIQUID WITH ALCOHOL ODOR

BOILING POINT: 181° F

FLASH POINT: 74° F

VAPOR DENSITY N/A

EVAPORATION RATE: N/A

SOLUBILITY IN WATER: COMPLETE

SPECIFIC GRAVITY (H₂O = 1): 0.892

IV. FIRE & EXPLOSION HAZARD DATA

FLASH POINT(Method Used): 74° F TAG CLOSED CUP : LEL: N/A UEL:N/A

EXTINGUISHING MEDIA: WATER, FOAM

SPECIAL FIRE FIGHTING PROCEDURES: TREAT PRODUCT AS FLAMMABLE

UNUSUAL FIRE AND EXPLOSION HAZARDS: NONE

V. REACTIVITY DATA

STABILITY: STABLE

CONDITIONS TO AVOID : NONE

INCOMPATIBILITY: ANIONIC SURFACTANT, SOAPS

HAZARDOUS DECOMPOSITION OR BYPRODUCTS: NONE KNOWN

POLYMERIZATION: WILL NOT OCCUR

CONDITIONS TO AVOID: NONE KNOWN

PRODUCT DOES NOT REACT WITH FORMALDEHYDE

MATERIAL SAFETY DATA SHEET

SUPER SANI CLOTH

VI. HEALTH HAZARD DATA

EFFECTS OF OVER EXPOSURE:

SKIN : MAY CAUSE SKIN IRRITATION
EYES : MAY CAUSE IRRITATION
INHALATION: NOT NORMAL ROUTE OF ENTRY
INGESTION: NOT NORMAL ROUTE OF ENTRY

EMERGENCY AND FIRST AID PROCEDURES:

SKIN CONTACT: IF RASH OR IRRITATION DEVELOPS, DISCONTINUE USE

EYE CONTACT: FLUSH WITH COOL WATER FOR 15 MINUTES

INHALATION: REMOVE TO FRESH AIR

INGESTION: CONSULT PHYSICIAN

TARGET ORGANS: NONE KNOWN

KNOWN CARCINOGEN: NONE KNOWN

VII. SPILL AND DISPOSAL PROCEDURE

SPILL CONTROL: PICK UP WITH ABSORBENT WASTE

DISPOSAL METHOD: TREAT AS FLAMMABLE LIQUID. OBEY LOCAL, STATE AND FEDERAL REGULATIONS

HANDLING AND STORAGE: AVOID HIGH HEAT AND SOURCES OF IGNITION

VIII. CONTROL MEASURES/PROTECTION

RESPIRATION: USE NIOSH APPROVED RESPIRATOR WHEN HANDLING BULK LIQUID IF EXPOSURE LEVEL IS EXCEEDED

VENTILATION: RECOMMENDED

PROTECTIVE GLOVES: RECOMMENDED

EYE PROTECTION: IF SPLASHING MAY OCCUR

HYGIENIC PRACTICES: FOLLOW GOOD HOUSEKEEPING PRACTICES

OTHER: NONE

IX. TRANSPORT/SHIPPING

DOT SHIPPING NAME: CONSUMER COMMODITY

TECHNICAL SHIPPING NAME: N/A

DOT SHIPPING CLASSIFICATION: ORM-D

DOT ID NO.: N/A

DOT LABEL REQUIREMENTS:

UN/NA NUMBER REGULATIONS: N/A

REPORTABLE QUANTITY: N/A

X. DISCLAIMER

THE INFORMATION FURNISHED HEREIN IS BELIEVED TO BE ACCURATE AND REPRESENTS THE BEST DATA CURRENTLY AVAILABLE TO US. NO WARRANTY, EXPRESSED OR IMPLIED IS MADE AND NICE-PAK PRODUCTS, INC. ASSUMES NO LEGAL RESPONSIBILITY OR LIABILITY RESULTING FROM ITS USE. ALL COMPONENTS OF THIS PRODUCT ARE LISTED UNDER T.S.C.A. INVENTORY.